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
FLEXIBLE SOLUTION

FOR PROXIMAL BONE PRESERVATION

The Epoch Hip Prosthesis is designed to address concerns about proximal femoral bone resorption resulting from stress shielding, as well as concerns about thigh pain that can result from a modulus mismatch between the prosthesis and natural bone, particularly with larger stems.

The design goal was to create an extensively porous-coated implant in sizes greater than 13mm without a significant increase in stem stiffness when compared to smaller sizes or proximally coated implants.

The larger Epoch Hip sizes (14-18mm) consist of an innovative composite construction that combines a tapered cobalt-chromium-molybdenum alloy core with an inert polyaryletherketone (PAEK) polymer matrix and a titanium fiber metal fixation surface. The bending stiffness of this low-modulus composite stem is 75% less than a

comparably shaped cobalt-chromium stem, and 50% less than a titanium alloy stem.¹ It more closely mimics the bending stiffness of bone, while maintaining performance characteristics comparable to conventional metal hip stems.¹

Studies suggest that smaller diameter stems, especially if manufactured from titanium alloy, do not produce a significant mismatch in flexural rigidity when compared to the natural femur.² The 11-13mm stems in the Epoch Hip family are manufactured from Titanium[®] Ti-6Al-4V Alloy, incorporating a circumferential proximal fiber metal pad and an adjacent corundumized surface.

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